Elementanalyse

TNO

ORGANISCH CHEMISCH INSTITUUT

juli 1975

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INSTITUTE FOR ORGANIC CHEMISTRY

july 1975

Organization for
Applied Scientific Research
in the Netherlands

Institute for Organic Chemistry

P.O. Box 5009
79 Croesestraat
Utrecht
the Netherlands
Telephone: 030-882721-extension 53

In charge of the
Element-analytical Section
W. J. Buis

Introduction

Because of an increasing interest in determinations of elements in organometallic and metal-organic compounds, modifications of the procedures for determination of C, H, Br, Cl, F, I, N, O and S have been worked out and possibilities to determine metals e.g. Cd, Co, Cu, Fe, Mg, Mn, Ni, Pb, Sb, Sn and Zn have been extended. A survey of the elements which can be determined is given on page 5.

Prices for single determinations are listed on page 6.

To ensure reliable results at least some knowledge of the composition and structure of the substance is essential. This knowledge determines the choice of procedure. Also important are the pre-treatment of the sample, and the proper choice of sample tube. Full instructions are given on page 4.

General instructions

Pretreatment of the sample

Often, the sample has to be dried before analysis. It is advisable to specify under what circumstances the sample has to be dried (e.g. P_2O_5 ; at 20° ; at 0.1 mm).

If no conditions for drying are stated, all samples will be analysed as received.

Bear in mind that large crystals occasionally include traces of the solvent. Due attention should also be given to the sample being homogeneous. Liquids and unstable or hygroscopic solids should preferably be mailed in relatively long, sealed tubes.

Stable and non-hygroscopic solids should be packed in carefully closed tubes.

Wrap corks in cigarette paper (against cork particles).

Provide each tube with a label in such a way that the sample contained in it is fully visible from the outside.

Information regarding the sample

Please, specify for each substance, if possible:

all elements present, this is necessary for the proper choice of destruction procedure;

melting or boiling point;

the approximate percentage of the elements to be determined (e.g. C 70–75%; H 7–9%; O 10–15%);

any special properties of the substance (e.g. hygroscopic, volatile, unstable, toxic, etc.);

whether treatment under nitrogen is necessary;

code number of the sample (also on the label).

Size of the sample

The quantity required for one single metal-determination is approximately 50 mg. For the determination of other elements 3–5 mg are required. Unless return is requested the remaining material will be destroyed after half a year.

Reference in case of publication

The following text is suggested to authors of papers who want to include an acknowledgement: Element analyses were carried out by the Element Analytical Section of the Institute for Organic Chemistry TNO, Utrecht, The Netherlands, under supervision of W. J. Buis.

SURVEY OF THE ELEMENTS WHICH CAN BE DETERMINED

The elements printed in red can be determined according to standard methods. For determination of one of the other elements, please apply to the Element-analytical Department.

H																	He
Li	Be											B	C	N	O	F	Ne
Na	Mg											Al	Si	P	S	Cl	Ar
K	Ca	Sc	Ti	V	Cr	Mn	Fe	Co	Ni	Cu	Zn	Ga	Ge	As	Se	Br	Kr
Rb	Sr	Y	Zr	Nb	Mo	Tc	Ru	Rh	Pd	Ag	Cd	In	Sn	Sb	Te	I	Xe
Cs	Ba	*	Hf	Ta	W	Re	Os	Ir	Pt	Au	Hg	Tl	Pb	Bi	Po	At	Rn
Fr	Ra	*	Ku														

Price-list

Effective from July 1st, 1975. Prices for a single determination are given in Dutch guilders. For a duplicate determination a discount of D fl 5.00 is calculated.

Carbon-Hydrogen	32.50
Bromine	29.00
Chlorine	29.00
Fluorine	45.00
Iodine	32.50
Nitrogen	27.00
Oxygen	31.00
Oxygen (addition of Induline base or "Merz" mixture) surcharge	10.00
Phosphorus	40.00
Sulphur	29.00
Boron	45.00
Mercury	75.00
Metal (as oxide or sulfate)	29.00
Metal (as residue with C-H determination)	17.50
Metal* (other methods)	55.00
Methoxy, Ethoxy	65.00
Handling in inert atmosphere (nitrogen)	10.00
Drying to constant weight	17.50

* Most elements can be determined by standard methods: Aluminium, Antimony, Arsenic, Barium, Bismuth, Cadmium, Calcium, Cobalt, Copper, Germanium, Gold, Indium, Iron, Lead, Magnesium, Manganese, Nickel, Silicon, Thallium, Tin, Titanium, Zinc. For substances containing elements which interfere upon analysis according to standard methods, supplementary charges will be made.

Please apply for discounts on contract or quantum basis.